


# Protocol: High Resolution Fixed Large Brain

## Purpose

High resolution, T2 weighted imaging of *ex vivo* larger brains (up to 8 cm long axis).



Large Brain Protocol – “C_FixedLargeBrainHighResolution”				
Type	Scan Name	In-plane resolution, slice thickness	Fov	Scan time
Localizer	1_tripilot			
T2 weighted	2_TurboRARE3D_117x469x625_1h	117 $\mu$ m x 469 $\mu$ m, 625 $\mu$ m	6cm x 6cm x 8cm	1 h
T2 weighted	3_TurboRARE3D_234x234x400_3.4h	243 $\mu$ m x 234 $\mu$ m, 400 $\mu$ m	6cm x 6cm x 8cm	2.5 h
T2 weighted	4_TurboRARE3D_117x234x305_3.4h	117 $\mu$ m x 234 $\mu$ m, 305 $\mu$ m	6cm x 6cm x 8cm	3.4 h

## Instructions

1. **Prepare the sample:** Place the sample on a sample holder that fits Coil A.
2. **Center the sample:** Center the sample in the RF coil center, then position the coil at isocenter in the magnet.
3. **Match, tune and set up scan as usual.**
4. **Load protocol:** Scans in this protocol are located in "C\_LargeBrain." After the localizer scan run the TurboRARE3D scan with the appropriate resolution as indicated in the table above.

## Notes

- Coil A, 72mm ID quad coil.

## Post-Processing Support

File conversion to DICOM, afni brik and niftii format.