

Bioengineering 280A
Principles of Biomedical Imaging

Fall Quarter 2009
MRI Lecture 2

TT Liu, BE280A, UCSD Fall 2009

Relaxation

An excitation pulse rotates the magnetization vector away from its equilibrium state (purely longitudinal). The resulting vector has both longitudinal M_z and transverse M_{xy} components.

Due to thermal interactions, the magnetization will return to its equilibrium state with characteristic time constants.

T_1 spin-lattice time constant, return to equilibrium of M_z

T_2 spin-spin time constant, return to equilibrium of M_{xy}

TT Liu, BE280A, UCSD Fall 2009


Longitudinal Relaxation

$$\frac{dM_z}{dt} = -\frac{M_z - M_0}{T_1}$$

After a 90 degree pulse $M_z(t) = M_0(1 - e^{-t/T_1})$


Due to exchange of energy between nuclei and the lattice (thermal vibrations). Process continues until thermal equilibrium as determined by Boltzmann statistics is obtained.

The energy ΔE required for transitions between down to up spins, increases with field strength, so that T_1 increases with \mathbf{B} .



TT Liu, BE280A, UCSD Fall 2009

T1 Values




Image, caption: Nishimura, Fig. 4.2

TT Liu, BE280A, UCSD Fall 2009

Transverse Relaxation

$$\frac{dM_{xy}}{dt} = -\frac{M_{xy}}{T_2}$$




Each spin's local field is affected by the z-component of the field due to other spins. Thus, the Larmor frequency of each spin will be slightly different. This leads to a dephasing of the transverse magnetization, which is characterized by an exponential decay.

T_2 is largely independent of field. T_2 is short for low frequency fluctuations, such as those associated with slowly tumbling macromolecules.

TT Liu, BE280A, UCSD Fall 2009

T2 Relaxation

Free Induction Decay (FID)




After a 90 degree excitation

$$M_{xy}(t) = M_0 e^{-t/T_2}$$

TT Liu, BE280A, UCSD Fall 2009

T2 Relaxation



TT Liu, BE280A, UCSD Fall 2009

Credit: Larry Frank

T2 Values

Tissue	T_2 (ms)
gray matter	100
white matter	92
muscle	47
fat	85
kidney	58
liver	43
CSF	4000


Table: adapted from Nishimura, Table 4.2

Solids exhibit very short T_2 relaxation times because there are many low frequency interactions between the immobile spins.

On the other hand, liquids show relatively long T_2 values, because the spins are highly mobile and net fields average out.

TT Liu, BE280A, UCSD Fall 2009

Example



T₁-weighted


Density-weighted

T₂-weighted

Questions: How can one achieve T₂ weighting? What are the relative T₂'s of the various tissues?

TT Liu, BE280A, UCSD Fall 2009

Example



(a) Four images, all obtained with a common TR=5 seconds and TE=90, 50, 20, 15 ms (shown in reading order).

(b) Six images obtained with a common TE=15 ms and TR=500, 1000, 2000, 3000, 4000, 5000 ms (shown in reading order).

Figure 8: Phantom data which illustrates signal intensity and contrast for bottles filled with jello of varying consistency. Where is T₁ long/short? How long, how short? The same for T₂? Which bottles might be pure water? Which jello is most firm? What pictures are the most T₁-, T₂- and PD-weighted?

TT Liu, BE280A, UCSD Fall 2009

Hanson 2009

Bloch Equation


$$\frac{d\mathbf{M}}{dt} = \underbrace{\mathbf{M} \times \gamma \mathbf{B}}_{\text{Precession}} - \underbrace{\frac{M_x \mathbf{i} + M_y \mathbf{j}}{T_2}}_{\text{Transverse Relaxation}} - \underbrace{\frac{(M_z - M_0) \mathbf{k}}{T_1}}_{\text{Longitudinal Relaxation}}$$

$\mathbf{i}, \mathbf{j}, \mathbf{k}$ are unit vectors in the x,y,z directions.

TT Liu, BE280A, UCSD Fall 2009

Free precession about static field

$$\begin{aligned} \frac{d\mathbf{M}}{dt} &= \mathbf{M} \times \gamma \mathbf{B} \\ &= \gamma \begin{vmatrix} \hat{i} & \hat{j} & \hat{k} \\ M_x & M_y & M_z \\ B_x & B_y & B_z \end{vmatrix} \\ &= \gamma \begin{pmatrix} \hat{i}(B_z M_y - B_y M_z) \\ -\hat{j}(B_z M_x - B_x M_z) \\ \hat{k}(B_y M_x - B_x M_y) \end{pmatrix} \end{aligned}$$



TT Liu, BE280A, UCSD Fall 2009

Free precession about static field

$$\begin{aligned} \begin{bmatrix} dM_x/dt \\ dM_y/dt \\ dM_z/dt \end{bmatrix} &= \gamma \begin{bmatrix} B_z M_y - B_y M_z \\ B_x M_z - B_z M_x \\ B_y M_x - B_x M_y \end{bmatrix} \\ &= \gamma \begin{bmatrix} 0 & B_z & -B_y \\ -B_z & 0 & B_x \\ B_y & -B_x & 0 \end{bmatrix} \begin{bmatrix} M_x \\ M_y \\ M_z \end{bmatrix} \end{aligned}$$

TT Liu, BE280A, UCSD Fall 2009

Precession

$$\begin{bmatrix} dM_x/dt \\ dM_y/dt \\ dM_z/dt \end{bmatrix} = \gamma \begin{bmatrix} 0 & B_0 & 0 \\ -B_0 & 0 & 0 \\ 0 & 0 & 0 \end{bmatrix} \begin{bmatrix} M_x \\ M_y \\ M_z \end{bmatrix}$$


Useful to define $M \equiv M_x + jM_y$

$$\begin{aligned} dM/dt &= d/dt(M_x + jM_y) \\ &= -j\gamma B_0 M \end{aligned}$$

Solution is a time-varying phasor

$$M(t) = M(0)e^{-j\gamma B_0 t} = M(0)e^{-j\omega_0 t}$$

Question: which way does this rotate with time?



TT Liu, BE280A, UCSD Fall 2009

Matrix Form with $B=B_0$

$$\begin{bmatrix} dM_x/dt \\ dM_y/dt \\ dM_z/dt \end{bmatrix} = \begin{bmatrix} -1/T_2 & \gamma B_0 & 0 \\ -\gamma B_0 & 1/T_2 & 0 \\ 0 & 0 & -1/T_1 \end{bmatrix} \begin{bmatrix} M_x \\ M_y \\ M_z \end{bmatrix} + \begin{bmatrix} 0 \\ 0 \\ M_0/T_1 \end{bmatrix}$$

TT Liu, BE280A, UCSD Fall 2009

Z-component solution

$$M_z(t) = M_0 + (M_z(0) - M_0)e^{-t/T_1}$$

Saturation Recovery

$$\text{If } M_z(0) = 0 \text{ then } M_z(t) = M_0(1 - e^{-t/T_1})$$

Inversion Recovery

$$\text{If } M_z(0) = -M_0 \text{ then } M_z(t) = M_0(1 - 2e^{-t/T_1})$$


TT Liu, BE280A, UCSD Fall 2009

Transverse Component

$$M \equiv M_x + jM_y$$

$$\begin{aligned} dM/dt &= d/dt(M_x + jM_y) \\ &= -j(\omega_0 + 1/T_2)M \end{aligned}$$


$$M(t) = M(0)e^{-j\omega_0 t} e^{-t/T_2}$$



TT Liu, BE280A, UCSD Fall 2009

Summary

- 1) Longitudinal component recovers exponentially.
- 2) Transverse component precesses and decays exponentially.




Source: <http://mrsrl.stanford.edu/~brian/mri-movies/>

TT Liu, BE280A, UCSD Fall 2009

Summary

- 1) Longitudinal component recovers exponentially.
- 2) Transverse component precesses and decays exponentially.



Fact: Can show that $T_2 < T_1$ in order for $|M(t)| \leq M_0$
Physically, the mechanisms that give rise to T_1 relaxation also contribute to transverse T_2 relaxation.


TT Liu, BE280A, UCSD Fall 2009

Gradients

Spins precess at the Larmor frequency, which is proportional to the local magnetic field. In a constant magnetic field $B_z = B_0$, all the spins precess at the same frequency (ignoring chemical shift).

Gradient coils are used to add a spatial variation to B_z such that $B_z(x, y, z) = B_0 + \Delta B_z(x, y, z)$. Thus, spins at different physical locations will precess at different frequencies.

TT Liu, BE280A, UCSD Fall 2009



Imaging: localizing the NMR signal

The local precession frequency can be changed in a position-dependent way by applying linear field gradients

$$\Delta B(x) \propto x$$

Resonant Frequency:
 $\nu(x) = \gamma B_0 + \gamma \Delta B(x)$

RF and Gradient Coils

Credit: R. Buxton

TT Liu, BE280A, UCSD Fall 2009

Gradient Fields

$$B_z(x, y, z) = B_0 + \frac{\partial B_z}{\partial x} x + \frac{\partial B_z}{\partial y} y + \frac{\partial B_z}{\partial z} z$$

$$= B_0 + G_x x + G_y y + G_z z$$

$$G_z = \frac{\partial B_z}{\partial z} > 0 \qquad G_y = \frac{\partial B_z}{\partial y} > 0$$

TT Liu, BE280A, UCSD Fall 2009

Interpretation

$$\Delta B_z(x) = G_x x$$

Spins Precess at $\gamma B_0 - \gamma G_x x$ (slower) Spins Precess at γB_0 Spins Precess at $\gamma B_0 + \gamma G_x x$ (faster)

TT Liu, BE280A, UCSD Fall 2009

Rotating Frame of Reference

Reference everything to the magnetic field at isocenter.



TT Liu, BE280A, UCSD Fall 2009


Spins



There is nothing that nuclear spins will not do for you, as long as you treat them as human beings.
Erwin Hahn


TT Liu, BE280A, UCSD Fall 2009

Phasors



TT Liu, BE280A, UCSD Fall 2009

Phasor Diagram



TT Liu, BE280A, UCSD Fall 2009

